# PUREZERO\* HG3 STERILE NITRILE GLOVES

12 in. / 30.5 cm

Hand Specific

Accelerator Free\*

**HALYARD\* PUREZERO\* HG3 Sterile Light Blue Nitrile Gloves**  
**HALYARD\* PUREZERO\* HG3 Sterile White Nitrile Gloves**

- 12 in./30.5 cm Length
- Accelerator Free \*
- Hand Specific Pairs
- Textured
- For the Sterile Critical Cleanroom Environment
- For Industrial Use Only

**NOTICE: THIS INSERT SHOULD BE FURNISHED OR MADE AVAILABLE TO THE USERS OF THESE GLOVES AS A SAFETY PRECAUTION.**

**This product is Category III PPE certified according to Regulation (EU) 2016/425 EEC.**

**Warning: Gloves offer protection against chemicals (Splash) and micro-organisms.**

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in the physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly.

For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Before use, inspect the gloves for defects or imperfections such as holes and tears. If the gloves are ripped or puncture during use, dispose of them immediately.

\* Not formulated with the following vulcanizing chemicals and accelerators: Sulfur, Thiurams, Thiazoles, Guanidines, Carbamates

**HALYARD\* PUREZERO\* HG3 Gants stériles en nitrile bleu clair**  
**HALYARD\* PUREZERO\* HG3 Gants stériles en nitrile blancs**

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www.halyardhealth.com/PPE-DOC-Info

O&M Halyard, Inc., 9120 Lockwood Blvd., Mechanicsville, VA 23116; In USA, please call 1-844-425-9273 • halyardhealth.com [CE] [RE] ArcRoyal Unlimited Company, Virginia Road, Kells; Co Meath, Ireland

O&M Halyard Belgium BVBA, Berkenlaan 8B, 1831 Machelen (Brab.), Belgium • Sponsored in Australia by O&M Halyard Australia Pty Ltd.; 52 Alfred Street S, Milsoms Point, NSW 2061

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Product certified to the PPE Regulation 2016/425 by Notified Body 2797: BSI Group The Netherlands B.V.; Say Building; John M. Keynesplein 9; 1066 EP Amsterdam

\*Registered Trademark or Trademark of O&M Halyard, Inc., or its affiliates. ©2018. All rights reserved. Patent and/or Pending Application. See www.halyardhealth.com/patents. 2020-05-07

15-HM-029-0-00