



**Cleanroom Gloves Product Information – HALYARD\* PUREZERO\* HG5 White Nitrile Gloves**

**Product Family: HALYARD\* PUREZERO\* HG5 White Nitrile Gloves**

**Product Codes: CLN30512XS, CLN30512SM, CLN30512MD, CLN30512LG, CLN30512XL**

	Permeation Test		Degradation Test
	EN 16523-1:2015+A1:2018		EN ISO 374-4:2019
Chemical	minimum breakthrough time (min)	Performance Level	Degradation (%)
Formaldehyde, 37%	>480	6	3.9
Hydrogen Peroxide, 30%	>480	6	5.4
Sodium Hydroxide, 40%	>480	6	7.1

AQL 1.0 Level 2 GI

Manufactured at our SAFESKIN\* Facility

Low Dermatitis Potential

**Product Family: HALYARD\* PUREZERO\* HG5 Blue SGX\* Nitrile Gloves**

**Product Codes: CLN92512XS, CLN92512SM, CLN92512MD, CLN92512LG, CLN92512XL**

	Permeation Test		Degradation Test
	EN 16523-1:2015+A1:2018		EN ISO 374-4:2019
Chemical	minimum breakthrough time (min)	Performance Level	Degradation (%)
Formaldehyde, 37%	>480	6	27.5
n-Heptane, 99%	>480	6	39.4
Sodium Hydroxide, 40%	>480	6	26.2

AQL 1.0 Level 2 GI

Manufactured at our SAFESKIN\* Facility

**Product Family: HALYARD\* PUREZERO\* HG5 Blue Nitrile 10 in. Gloves**

**Product Codes: CLN90511XS, CLN90511SM, CLN90511MD, CLN90511LG, CLN90511XL**

	Permeation Test		Degradation Test
	EN 16523-1:2015+A1:2018		EN ISO 374-4:2019
Chemical	minimum breakthrough time (min)	Performance Level	Degradation (%)
Formaldehyde, 37%	>480	6	57.1
Hydrogen Peroxide, 30%	>480	6	39.2
Sodium Hydroxide, 40%	>480	6	11.5

AQL 1.0 Level 2 GI

Manufactured at our SAFESKIN\* Facility