

**Informace o sterilizačním sáčku**


Sterilizační cykly pro 48495, 48521, and 48540, 48522, 48523, 48557, 48556, 48555, 48524, 48506, 48549, 48499, 48510, 48518, 48550, 48498, 48526, 48519, 48513, 48541, 48546, 48542, 48497, 48538, 48507, 48512, 48558, 48508, 48525, 48501, 48520, 48500, 48496, 48505, 48494, 48503



Pára nebo chemická pára na fialovou



Etylenoxidový plyn na hnědou



Sterilizační cyklus pro 48534, 48535, 48531, 48532



Plazmový proces s  
H<sub>2</sub>O<sub>2</sub> na růžovou

