

**Produktinformation vedrørende sterilisationsposer**


Steriliseringscyklusser til 48495, 48521, and 48540, 48522, 48523, 48557, 48556, 48555, 48524, 48506, 48549, 48499, 48510, 48518, 48550, 48498, 48526, 48519, 48513, 48541, 48546, 48542, 48497, 48538, 48507, 48512, 48558, 48508, 48525, 48501, 48520, 48500, 48496, 48505, 48494, 48503



Damp eller kemisk  
damp til lilla



E.O. Gas til brun



Steriliseringscyklus til for 48534, 48535, 48531, 48532



H<sub>2</sub>O<sub>2</sub> plasmaproces til lyserød

