

Sterilization Pouch Product Information


Sterilization Cycles for 48495, 48521, 48540, 48522, 48523, 48557, 48556, 48555, 48524, 48506, 48549, 48499, 48510, 48518, 48550, 48498, 48526, 48519, 48513, 48541, 48546, 48542, 48497, 48538, 48507, 48512, 48558, 48508, 48525, 48501, 48520, 48500, 48496, 48505, 48494, 48503



STEAM OR CHEMICAL VAPOR to PURPLE



E.O. GAS to BROWN



Sterilization Pouch Product Information

Sterilization Cycle for 48534, 48535, 48531, 48532



**H₂O₂ PLASMA PROCESS
to PINK**

