

Informazioni sulla busta di sterilizzazione


Cicli di sterilizzazione per 48495, 48521, 48540, 48522, 48523, 48557, 48556, 48555, 48524, 48506, 48549, 48499, 48510, 48518, 48550, 48498, 48526, 48519, 48513, 48541, 48546, 48542, 48497, 48538, 48507, 48512, 48558, 48508, 48525, 48501, 48520, 48500, 48496, 48505, 48494, 48503



Con il vapore o vapore chimico il colore cambia in porpora



Con il gas EO il colore cambia in marrone



Ciclo di sterilizzazione per 48534, 48535, 48531, 48532



Con il trattamento a plasma di
H₂O₂ il colore cambia in rosa

