

Informações sobre a bolsa de esterilização


Ciclos de esterilização para 48495, 48521, 48540, 48522, 48523, 48557, 48556, 48555, 48524, 48506, 48549, 48499, 48510, 48518, 48550, 48498, 48526, 48519, 48513, 48541, 48546, 48542, 48497, 48538, 48507, 48512, 48558, 48508, 48525, 48501, 48520, 48500, 48496, 48505, 48494, 48503



Vapor ou vapor químico para roxo



Gás OE para castanho



Informações sobre a bolsa de esterilização

Ciclo de esterilização para 48534, 48535, 48531, 48532



Processo de plasma H₂O₂ para rosa

