

灭菌袋产品信息


的 灭菌 周期 48495, 48521, 48540, 48522, 48523, 48557, 48556, 48555, 48524, 48506, 48549, 48499, 48510, 48518, 48550, 48498, 48526, 48519, 48513, 48541, 48546, 48542, 48497, 48538, 48507, 48512, 48558, 48508, 48525, 48501, 48520, 48500, 48496, 48505, 48494, 48503



蒸气或化学蒸气至紫色



E.O.气体至棕色



的 灭菌 周期 48534, 48535, 48531, 48532



H₂O₂ 等离子体处理至粉红色

